



## Teeth Whitening Pro Kit

Art no: 125051 & 125052

Up to 7 shades whiter teeth in 10 days

### Claims:

- Triple Boost Max Technology that cleans with detergents and whitens with sodium bicarbonate and PAP for increased efficiency.
- Safe and without side effects
- Kit includes Protect & Repair gel that protects against new stains and re-mineralizes enamel.
- Customized and reusable mouth trays
- Up to 10 days treatment

### Product description:

The Pro Kit is a complete system to whiten teeth at home fast, easy and gentle. It contains 10 ml Peroxide Free Triple Boost Max Gel, which is used together with the thermoform mouth trays that are individually fitted using hot water.

The Pro Kit also contains a Protect & Repair gel which is used post whitening for protection and strengthening of teeth.

### Working mechanism in short:

#### The Triple Boost Max gel

The gel contains PAP ( $\epsilon$ -phthalimido peroxy hexanoic acid) + sodium bicarbonate + coco sulfate which work together for the most effective peroxide free whitening treatment available. The PAP ingredient is an organic peracid resulting in a high degree of oxidation. Oxidation is necessary for the bleaching procedure as it neutralizes organic double bonds that cause dental discolorations. However, the activity of PAP is not based on the release of hydrogen peroxide and the gel is 100% peroxide free.

A fourth whitening effect can be created with the Beconfident<sup>®</sup> LED light which accelerates the result significantly using ferrous gluconate activator and blue light radiation.

#### The Protect and Repair Gel (10 ml tube)

The enamel consists of about 97% of hydroxyapatite, which is why we added Hydroxyapatite to support the repair process after whitening. Even incipient caries lesions can be repaired. Daily brushing with a hydroxyapatite-containing toothpaste also helps to seal the so-called dentine tubules. In the case of exposed tooth necks, the pain sensitivity associated with thermal and chemical stimuli can be significantly reduced.

The gel also contains a sodium salt called Tetranatrium pyrophosphat that helps to dissolve the discoloration of the teeth and to remove them more easily. Furthermore, the pyrophosphate has the property to restrict the formation of tartar. It also contains sorbitol which has moisturizing characteristics

Group Beconfident Europe AB

Tennisvägen 31, SE-26933 Båstad, Sweden

E-mail: [info@beconfident.info](mailto:info@beconfident.info) Phone: +46-(0)702 56 50 00 VAT no: SE 556853-0041-01

Corporate site: [www.beconfident.info](http://www.beconfident.info), [www.beconfident.co.uk](http://www.beconfident.co.uk), [www.beconfident.se](http://www.beconfident.se) & many more

Beconfident- Top Quality Teeth whitening and Skincare systems are designed in Sweden & Made in USA and Sweden.

Beconfident<sup>®</sup> + WhiteAmin<sup>®</sup> + SimpleSmile<sup>®</sup> are registered trademarks in Europe and Beconfident<sup>®</sup> is a registered trademark in United States of America. IOS and Android Intelligent Apps belong to Beconfident Skincare AB. Other brands are pending.





and Sodium Fluoride to prevent caries. The Chamomile is an extract from the flowers of chamomile which has a soothing effect on the gums. The gel contains Sodium monofluorophosphate (920 ppm F-).

#### Unique Competitive Advantages Triple Boost Max gel:

- Unique Triple Boost Max technology for a safe process and efficient result that uses 3 different mechanisms to whiten teeth
- Whitening treatments last longer because stains are actually removed.
- Individually fitted to optimize the contact between the gel and tooth surface.
- Non-Abrasive. Does not scratch teeth enamel and works safely for all users.
- Includes free Teeth Whitening app to guide and support customer during treatment in order to achieve the best result.

#### Unique Competitive Advantages Protect & Repair Gel:

- Protects teeth and prolongs results after whitening
- Repairs and strengthens damaged enamel
- No need for trays, applies easily in 30 seconds without any hassle
- Immediate optic whitening effect

#### Ingredients:

##### Triple Boost Teeth Whitening Gel:

Glycerin, Sodium Bicarbonate, Phthalimidoperoxycaproic acid, Mentha Piperita Oil, Sodium Coco-Sulfate, Xylitol, Chondrus Crispus Powder, Sorbitol, Aqua, Aloe Barbadensis Leaf Juice, Punica Granatum Extract, Chamomilla Recutita Flower Oil, Ferrous Gluconate, Limonene, Linalool.

##### Protect & Repair Teeth Gel:

Aqua, Glycerin, Sorbitol, Calcium Hydroxyapatite, Carbomer, Tetrasodium Pyrophosphate, Aloe Barbadensis Leaf Juice, Cocos Nucifera Oil, Sodium Coco-Sulfate, Petrolatum, Sodium Monofluorophosphate, Sodium Saccharin, Xanthan Gum, Chamomilla Recutita Flower Extract, Mentha Piperita Oil, Tocopheryl Acetate, Polysorbate 20, Phenoxyethanol, Triethanolamine, CI 42090.

#### The Pro kit contains the following items:

- 1) 1 x 10ml Syringe with Triple Boost Max Gel
- 2) 1 x 10 ml Protect & Repair
- 3) 2 shapeable mouth trays
- 4) Gel applicator pen
- 5) Instruction sheet and a shade guide.

Beconfident is not conducting any testing on animals.

Group Beconfident Europe AB

Tennisvägen 31, SE-26933 Båstad, Sweden

E-mail: [info@beconfident.info](mailto:info@beconfident.info) Phone: +46-(0)702 56 50 00 VAT no: SE 556853-0041-01

Corporate site: [www.beconfident.info](http://www.beconfident.info), [www.beconfident.co.uk](http://www.beconfident.co.uk), [www.beconfident.se](http://www.beconfident.se) & many more

Beconfident- Top Quality Teeth whitening and Skincare systems are designed in Sweden & Made in USA and Sweden.

Beconfident® + WhiteAmin® + SimpleSmile® are registered trademarks in Europe and Beconfident® is a registered trademark in United States of America. IOS and Android Intelligent Apps belong to Beconfident Skincare AB. Other brands are pending.



## Product Data:

Art No	Dimension HxWxD (mm)	Gross Weight	Net Weight	Units Display Box	Units mastercart.	Unit EAN
125051 & 125052	175x120x33	153g	121g	4	102	See below
CPNP Code	GMP ISO 22716	Country of Origin	Quality Test Pharmacies	Box bio-degradable	Pallet	RRP Euro incl VAT(+/-)
	Follows	USA		Yes	12	125 €uro

## Article numbers & EAN Codes

EU NORTH -51 (EN, SV, NO/DK, FI): 7350064167021
EU CENTRAL -52 (EN, FR, NL, DE, IT, ES): 7350064168011

## Group Beconfident Europe AB

Tennisvägen 31, SE-26933 Båstad, Sweden

E-mail: [info@beconfident.info](mailto:info@beconfident.info) Phone: +46-(0)702 56 50 00 VAT no: SE 556853-0041-01

Corporate site: [www.beconfident.info](http://www.beconfident.info), [www.beconfident.co.uk](http://www.beconfident.co.uk), [www.beconfident.se](http://www.beconfident.se) & many more

Beconfident- Top Quality Teeth whitening and Skincare systems are designed in Sweden & Made in USA and Sweden.

Beconfident® + WhiteAmin® + SimpleSmile® are registered trademarks in Europe and Beconfident® is a registered trademark in United States of America. IOS and Android Intelligent Apps belong to Beconfident Skincare AB. Other brands are pending.